

Redefine Group A Strep Test Performance at the Point of Care

Group A Strep (GAS) – A Common Bacterial Cause of Illness

- **Most common bacterial cause:** GAS is responsible for 5%–15% of sore throat visits in adults and 20%–30% in children¹
- **Clinical diagnosis alone is unreliable:** Signs and symptoms of GAS and non-streptococcal pharyngitis overlap so broadly that accurate diagnosis based on clinical grounds alone is usually impossible¹
- **Testing better enables antimicrobial stewardship:** As many as 10 million antibiotic prescriptions per year are directed toward respiratory conditions for which they are unlikely to provide benefits²

The BD Veritor™ System – The First CLIA-waived Digital Immunoassay (DIA) for the Rapid Detection of Group A Strep (GAS) With an Instrumented Result

- This digital, rather than visual, test result provides greater consistency regardless of the user's experience
- Reliable results available in minutes
- High sensitivity and specificity performance was established vs bacterial culture in a multicenter clinical trial (N=692)

Sample Type	Sensitivity*	Specificity*
Throat swab sample	95.4% (95% CI: 90.3%, 97.9%)	95.7% (95% CI: 93.7%, 97.1%)

*Reference method: bacterial culture; data from package insert

Group A
STREP



BD Veritor™ System

Changing the Way You View Rapid Testing

Redefine Group A Strep Test Performance at the Point of Care

Streamlined Workflow – Provides a digital result in minutes



Easy sample processing
Unitized tube containing the correct volume of process reagent facilitates workflow. Processing requires addition of 3 drops of Reagent 1 and 1-2 minutes incubation



Ready in minutes
Test device is ready to insert into reader 5 minutes after sample is added



Insert and read
Simple one-touch button readies the reader for test device insertion



Results delivered
Once the test device is inserted in the reader, an objective, digital test result is displayed in 10 seconds

Ordering Information

					
Description	Cat. No.	Qty.	Description	Cat. No.	Qty.
BD Veritor™ System Reader	256055	1	BD Veritor™ System Group A Strep CLIA-walved Kit	256040	30



BD Diagnostics
www.bd.com/ds

References: 1. Infectious Diseases Society of America. Clinical Practice Guidelines for the Diagnosis and Management of Group A Streptococcal Pharyngitis: 2012 Update by the Infectious Diseases Society of America. *Clin Infect Dis.* 2012. doi:10.1093/cid/cis629. 2. Hersh AL, et al. Principles of judicious antibiotic prescribing for bacterial upper respiratory tract infections in pediatrics. *Pediatrics.* doi:10.1542/peds.2013-3260.

BD, BD Logo and BD Veritor are trademarks of Becton, Dickinson and Company
©2014 BD 223679 February 2014

BD Veritor™ System

Changing the Way You View Rapid Testing