

WHEATON® CELLINE™ BIOREACTORS

Multiuse Membrane Culture Flasks for
Antibody and Protein Production



DURAN
WHEATON
KIMBLE


Excellence in your hands

 **fisher scientific**
part of Thermo Fisher Scientific

WHEATON®

CELLine™ Bioreactors

Advanced Cell Cultivation Devices



In an effort to push science forward, DWK Life Sciences is manufacturing a new line of membrane driven flasks for high density cell culture. The flasks are designed to enhance small scale bio-production for antibody and protein generation. Conventional in vivo or in vitro cell culture methods can be laborious, may have low cell density and require significant purification. CELLine flasks address these three areas of limitation observed in static tissue culture flasks.

Handling Requirements — CELLine flasks reduce the handling requirement by requiring less consumable items and allowing longer run times due to their unique metabolite regulating upper membrane. This membrane allows for bulk media storage during operation to ensure constant and regulated nutrient access for the cells.

Cell Densities — CELLine flasks ensure maximum gas exchange by placing the gas permeable lower membrane directly next to the cells. This allows for optimal oxygen and carbon dioxide transfer for metabolizing cells.

Purification — The upper and lower membranes form an optimized compartment for cell proliferation. This allows for the reduced use of growth factors and hormones and concentrates the antibody and proteins of interest.

Benefits of CELLine Flask

- Multiuse flask
- High cell density and high product concentration
- Reduces operation time
- Decreases use of consumables
- Cost-efficient, space saving, and stackable
- No additional equipment required for operation


How does the CELLine Flask Work?

Media Compartment — The media compartment allows for bulk storage of cell culture growth medium. This reduces the media refreshing requirement significantly as the media compartment is fifty times the size of the cell compartment.

Metabolite Regulating Upper Membrane — The upper dialysis membrane has a 10 kDa cut off limit. This regulates the flow of metabolites to and from the cell compartment and retains all proteins in the cell compartment.

Cell Compartment — The cell compartment provides the ideal area to inoculate and achieve high density cultures. The compartment concentrates cells, their products, and limits the requirement for any exogenous growth factors.

Gas Permeable Lower Membrane — With static cultures, gas transfer rates can be the limiting factor in high density cultures. The CELLine flask places the cells directly against the gas permeable membrane to achieve optimal levels of oxygen and carbon dioxide.



Exploded view of device; unit is packaged fully assembled

WHEATON®


CELLine™ Bioreactors Advanced Cell Cultivation Devices



1. Warm nutrient medium and prepare upper membrane by adding a small amount of media to the media compartment.



2. Inoculate the cell compartment of the device.



3. Fill the media compartment.



4. Incubate in a CO₂ incubator for approximately 3-7 days depending on the optimum harvesting schedule.




5. Change media from the media compartment and/or harvest the cells from the cell compartment.

- a. The higher the cell density, the more protein production. Higher cell densities will require more frequent media changes.

- b. Recommended 3, 5 and 7 day feeding schedules are available. Schedules can be tailored to meet your production and timing requirements.



6. Incubate flask and repeat the media change/harvesting cycle until desired amount of protein is generated.



WHEATON®

CELLine™ Bioreactors

Advanced Cell Cultivation Devices

Ordering Information

Cat. No.	Flask Type	Culture Type	Media (mL) Compartment Size	Cell (mL) Compartment Size	Qty/ Case
02-912-744	CELLine 1000	Suspension	1000	15	1
02-912-741	CELLine 1000	Suspension	1000	15	3
02-912-745	CELLine 1000-AD	Adherent	1000	15	1
02-912-742	CELLine 1000-AD	Adherent	1000	15	3
02-912-746	CELLine 350	Suspension	350	5	1
02-912-743	CELLine 350	Suspension	350	5	5

Additional Information

Seeding And Harvesting Densities	CELLine 350	CELLine 1000
Preculture (Viable Cells)	7.5 x 10 ⁶	22.5 x 10 ⁶
Inoculation Volume (mL)	5	15
Inoculation Concentration (Viable Cells/mL)	1.5 x 10 ⁶	1.5 x 10 ⁶
Harvesting Concentration (Viable Cells/mL)	20-40 x 10 ⁶	20-40 x 10 ⁶
Titer (mg/mL)	1 - 10	1 - 10
Antibody Yield per Month (mg)	20-200	60-600



WHEATON® and CELLine™ are trademarks or registered trademarks of DWK Life Sciences LLC. Copyright © 2019 DWK Life Sciences LLC. Printed in USA.
U.S. Patent: 5,693,537; Canada Patent: 2,193,810; Europe Patent: 0 769 048; Japan Patent: 3608664

FBRO_003_rv2 05/2018

In the United States:

For customer service, call 1-800-766-7000

To fax an order, use 1-800-926-1166

To order online: fishersci.com

In Canada:

For customer service, call 1-800-234-7437

To fax an order, use 1-800-463-2996

To order online: fishersci.ca



© 2018 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at fishersci.com/trademarks.