

## SAFETY PACK LABGLASS


[www.fisheredu.com/eisco](http://www.fisheredu.com/eisco)

- Convenient sets of high-quality borosilicate 3.3 glassware
- Packaged in thick, custom-fit foam inserts for safe transport and storage
- Superior heat tolerance and chemical resistance, autoclavable at 121°C
- Volumetric glassware is hand calibrated in Eisco's ISO 17025 compliant laboratory



### Eisco™ Safety Pack Mixed Glassware Set - Includes:

- 50ml, 100ml & 250ml Beakers with white enamel graduations
- 10ml, 50ml & 100ml Class-A, ASTM Graduated Cylinders featuring hexagonal bases, plastic protective collars, and blue enamel graduations
- 50ml, 100ml & 250ml Conical Flasks with white enamel graduations



Eisco™ Safety Pack Conical Flask Sets	
Code	Includes
S81267	50ml, 100ml, 250ml
S81268	100ml, 250ml, 500ml
S81269	500ml, 1000ml, 2000ml
S81270	50ml, 100ml, 250ml, 500ml, 1000ml



Eisco™ Safety Pack Volumetric Flask Sets	
Code	Includes
S81271	10ml, 25ml, 50ml Class-A, ASTM
S81272	50ml, 100ml, 250ml Class-A, ASTM



Eisco™ Safety Pack Graduated Cylinder Sets	
Code	Includes
S81273	10ml, 25ml, 50ml Class-A, ASTM
S81274	25ml, 50ml, 100ml Class-A, ASTM
S81275	50ml, 100ml, 250ml Class-A, ASTM



Eisco™ Safety Pack Beaker Sets	
Code	Includes
S27113	50ml, 100ml, 250ml

Distributed by Fisher Science Education. Contact us today:

**In the United States**

Order online: [fisheredu.com](http://fisheredu.com)

Call customer service: 1-800-955-1177

**Fisher Science Education**

4500 Turnberry Drive

Hanover Park, IL 60133

 **fisher science**  
education  
part of Thermo Fisher Scientific