



Fisherbrand Traceable™: Setting the Standard for Temperature Monitoring in Clinical Trials

Simple and accurate digital data logging and remote monitoring

USB Data Retrieval

Download data from device to USB flash drive

Wi-Fi/Ethernet Connected

Real-time, remote monitoring—receive alerts on mobile device or email



USB data logging

- Log data continuously
- Download memory data to PC

Cloud service monitoring

- 24/7 remote monitoring
- Built-in VFC reporting

	Traceable™ Excursion-Trac™ USB Data Logging	Traceable™ Wi-Fi or Ethernet Data Logger (with TraceableLIVE™ Cloud Service)
DATA STORAGE	2 years (resettable)	Unlimited cloud storage
USER INTERFACE	Local display	Local display, Web, smartphone
ALARMS	Audible, visual	Email, SMS text, Push notification, audible, visual
DATA OUTPUT	USB data stick to Excel file	Real-time graphs, reports, alarm history, statistics, audit trail, cal due notice
EMAIL/TEXT ALERTS	✗	✓
NIST-TRACEABLE CALIBRATION CERT	✓	✓
VACCINE FOR CHILDREN COMPLIANT	✓	✓
21 CFR 11 COMPLIANT OPTIONS	✗ (available on Memory-Loc options -see fishersci.com)	✓
APPLICATION	Cat. No.	Cat. No.
FRIDGE/FREEZER	15-081-123	Wi-Fi: 15-079-615 Ethernet: 15-079-704
ULTRA-LOW FREEZER -80°C	15-079-624	Wi-Fi: 15-079-626 Ethernet: 15-079-705
LIQUID NITROGEN DEWAR	15-079-669	15-079-639
ROOM MONITORING & INCUBATORS		
TEMPERATURE	15-081-123	15-079-616
TEMPERATURE/HUMIDITY	15-079-665	15-079-628
TEMPERATURE/HUMIDITY/PRESSURE	15-079-666	15-079-664
CO₂/INCUBATOR	✗	15-079-662

In the United States

Order online: fishersci.com
 Fax an order: 1-800-926-1166
 Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca
 Fax an order: 1-800-463-2996
 Call customer service: 1-800-234-7437

