




Innovative Testing Solutions for Your Most Pressing Concerns

- Help ensure quality and safety
- Find rapid and reliable testing essentials
- Meet regulatory requirements
- Avoid supply chain gaps with choice and convenience

The North American functional beverage market is expected to exceed **\$44 million** by 2025 and will require additional life sciences-based quality control and safety testing to keep consumers safe.¹



North America leads the global food and beverage market and is expected to grow to over **\$93 million** at a CAGR of

7.8% from 2022 to 2025²

By **2025**, the global beverage additives market is expected to be worth

\$11 million³

Confidence at Every Step

Beverage Quality Testing Solutions to Satisfy Evolving Tastes

Access over **1.7 million** food and beverage testing products, including chemicals, equipment, consumables, and more, from close to **6,500** suppliers.

The CBD-infused beverage market is evolving and diversifying and is expected to reach **\$1.7 billion** in the U.S. by 2025.⁴



Confidence at Every Step

Beverage Quality Testing Solutions
to Satisfy Evolving Tastes



**Millipore
Sigma**

MilliporeSigma EZ-Fit Filtration Products

Reduce the risk of sample contamination with MilliporeSigma™ EZ-Fit™ products. Their performance, ease of use, and sustainability support your beverage microbiological testing.



Choose 1, 3, or 6 Heads

EZ-Fit Manifolds

Remove EZ-Fit manifold filtration heads, valves, and tubing connections without additional tools. The low height makes them more suitable for use inside laminar flow hoods. Filtration heads for specific filtration devices sold separately.



Quiet, Efficient Operation

EZ-Stream Vacuum Pump


Save space and reduce liquid waste handling with the compact MilliporeSigma™ EZ-Stream™ Vacuum Pump, which produces pressures up to 600 to 700 mbar, is quiet during use (60 dB), and no water trap is needed.



Improved Microbial Recovery

EZ-Fit Funnels

Reliably recover microbes in bioburden testing with disposable pre-sterilized EZ-Fit filtration units. Simply collapse the funnel and add liquid media to use like an agar plate or use the protective lip to easily remove solid media.



Designed for Sustainability

EZ-Fit Manifolds

Help meet your sustainability goals with EZ-Fit manifolds, which are designed to help reduce environmental and health impacts, energy consumption, and waste.



Browse our full collection of beverage testing products at fishersci.com/food-beverage-testing and fishersci.ca/food-beverage-testing.

1. Functional Foods and Beverages: Global Markets, BBC Research (2020, December), page 94. <https://www.bccresearch.com/market-research/food-and-beverage/functional-food-market.html>
2. Functional Foods and Beverages: Global Markets, BBC Research (2020, December), page 86. <https://www.bccresearch.com/market-research/food-and-beverage/functional-food-market.html>
3. Global Market for Food Additives, BBC Research (2021, May), page 22. <https://www.bccresearch.com/market-research/food-and-beverage/food-additives-market.html>
4. Cannabis: CBD-infused Beverages: Looking at the sector's next frontier, Credit Suisse Equity Research (2019, December 5).

Distributed by Fisher Scientific. Contact us today:

In the United States

Order online: fishersci.com
Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca
Call customer service: 1-800-234-7437

 **fisher scientific**
part of Thermo Fisher Scientific