

# Low-Form Polypropylene Beakers

Tapered design  
facilitates easy  
handling



# Low-Form Beakers

Autoclavable and shatter resistant



Smoothly tapered Fisherbrand™ beakers made from virgin, high clarity polypropylene. This extensive range of plastic beakers from 25ml to 5 liters is break resistant and durable.

- Resistant to most acids, bases, and many commonly used solvents
- Anti-lock feature prevents beakers from sticking together when stacked
- Withstand temperatures up to 121°C (250°F)
- Easy-to-read, blue printed graduations
- “No Flame” symbol and SPI recycling code

## Ordering Information

Product	Cat. No.
<b>Fisherbrand Low-Form Polypropylene Beakers</b>	
25mL PP Tapered Beaker, 32mm x 51mm h 10/Pack	02-591-25
50mL PP Tapered Beaker, 39mm x 61mm h 10/Pack	02-591-26
100mL PP Tapered Beaker, 48.5mm x 73mm h 10/Pack	02-591-27
250mL PP Tapered Beaker, 65.5mm x 96mm h 10/Pack	02-591-28
400mL PP Tapered Beaker, 74.5mm x 111mm h 10/Pack	02-591-29
500mL PP Tapered Beaker, 84mm x 121mm h 10/Pack	02-591-30
600mL PP Tapered Beaker, 84mm x 126mm h 5/Pack	02-591-31
1000mL PP Tapered Beaker, 99mm x 148mm h 5/Pack	02-591-32
2000mL PP Tapered Beaker, 125mm x 188mm h 1/Pack	02-591-33
3000mL PP Tapered Beaker, 143mm x 214mm h 1/Pack	02-591-34
5000mL PP Tapered Beaker, 164mm x 250mm h 1/Pack	02-591-35

Find out more at [fishersci.com](https://fishersci.com)

### In the United States:

For customer service, call 1-800-766-7000

To fax an order, use 1-800-926-1166

To order online: [fishersci.com](https://fishersci.com)

### In Canada:

For customer service, call 1-800-234-7437

To fax an order, use 1-800-463-2996

To order online: [fishersci.ca](https://fishersci.ca)

