

Protect

Nalgene face shield

Thermo Scientific™ Nalgene™ Polycarbonate Face Shields provide

- Optically clear shield
- Protection from UV light (200–360 nm range)
- Biohazard and transparent
- Complies with OSHA 29 CFR 1910.1030 and ANSI 287.1-1989 standards for use with blood-borne pathogens

Ordering information

Description	No. per case	Cat. No.
Nalgene Polycarbonate Face Shield, adjustable, 11.6 x 7.8 in. (l x w), clear	4	11-409-5

Do not use with radioisotopes. Use Thermo Scientific™ Nalgene™ beta-radiation shields (Cat. No. 6700-1812 or 6700-2418) when working with beta-radioisotopes. Do not autoclave.

This product is intended for General Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customers' specific use or application.

Find out more at fishersci.com

Contact us today:

In the United States

Order online: fishersci.com
 Fax an order: 1-800-926-1166
 Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca
 Fax an order: 1-800-463-2996
 Call customer service: 1-800-234-7437

