

Professional Installations

from the Service Partner You Trust


You can't go wrong with a professional installation package from METTLER TOLEDO.

Our IPacs provide you with a professionally documented installation and an operational qualification to ensure your equipment will function according to its operational specifications in your environment.

StarterPac Installation

Start working right away and be sure your equipment is properly configured with the StarterPac solution. This package is a low cost option that can be used in any industry and includes:


- Documented installation by a factory-trained technician
- Efficient start-up procedures
- Certificate of NIST traceability
- Basic user training



IPac Standard Qualification

The IPac Standard Qualification is the best choice to ensure accuracy and meet national and international Quality Management requirements such as the ISO 9000 series and food standards. This installation package applies to all industries and includes:


- Professional installation by a factory-trained technician
- Initial installation/qualification documentation (IQ/OQ)
- Easy integration into existing quality system
- Basic user training
- Verification of performance tailored to your specific quality requirements
- Choice of ISO17025 calibration report or ISO900x manufacturers calibration certificate



Moisture IPac Standard Qualification

The Moisture IPac is a documented comprehensive installation and initial calibration of your Moisture Analyzer. This Pac includes:

- Professional installation by a factory-trained technician
- Initial installation/qualification documentation (IQ/OQ)
- Easy integration into existing quality system
- Basic user training
- Verification of performance tailored to your specific quality requirements
- Temperature calibration
- ISO900x manufacturers calibration certificate



METTLER TOLEDO

Equipment Qualifications

From the Service Partner You Trust

Ideal for customers in the pharmaceutical, biotechnology, and food industries, METTLER TOLEDO's equipment qualification packages are perfect for customers who prefer the entire qualification process be documented in a format designed and provided by the manufacturer.

EQPac Comprehensive Qualification Includes


IQ (Installation Qualification): Includes Installation and Adjustment of leveling and filter settings

OQ (Operational Qualification): Includes Testing according to SOP, Short Introduction to equipment per operating instructions and Calibration Certificate

PQ (Performance Qualification): Includes Support with performance qualification and overall documentation in METTLER TOLEDO Log Book.

Our EQ Pacs for both Precision and Analytical balances include the following benefits:

- Professional installation by a factory-trained technician
- Initial installation/qualification documentation (IQ/OQ)
- Documented OQ, PQ and MQ process
- Easy integration into existing quality system
- Full testing scope
- Basic user training
- Verification of performance tailored to your specific quality requirements
- Choice of ISO17025 calibration report, ISO900x manufacturers calibration certificate or measurement uncertainty USP certificate



Equipment Qualification Available

- Analytical and micro balances
- Precision balances
- Moisture analyzers*

*Only certificate option available: ISO900x manufacturers calibration certificate

METTLER TOLEDO